

DupLEX-ATM
YEAST
TWO-HYBRID

SYSTEM

User's Manual
Version 2.7
3/98

OriGene Technologies, Inc.
13 Taft Court, Suite 111
Rockville, MD 20850
1-888-267-4436
FAX: 301-340-9254

DupLEX-A Yeast Two-Hybrid System

Index

<u>I.</u>	Introduction	3
<u>II.</u>	List of components, additional materials required	4-7
	A. Yeast strains	4
	B. Reporter gene (lacZ) plasmids	4
	C. Bait plasmids	4
	D. Target plasmid	5
	E. Control plasmids	5
	F. Primers	5
	G. Other items included in the kit	6
	H. Materials required that are not included in the kit	6-7
<u>III.</u>	Recipes for media and solutions	8-11
	A. Yeast growth media	8-9
	B. Bacterial growth media	9-10
	C. Solutions	10-11
<u>IV.</u>	Library information	12
<u>V.</u>	Working with yeast	12
<u>VI.</u>	Flowchart of a complete library screen	13
<u>VII.</u>	Screening libraries	14-22
	A. Constructing the bait	14
	B. Testing the autoactivation potential of the bait (lacZ reporter)	14-15
	C. Testing the autoactivation potential of the bait (LEU2 reporter)	15-16
	D. Testing the bait's ability to enter the nucleus and bind LexA operators	16-17
	E. Performing a large-scale library transformation	17-18
	F. Screening for potential positive transformants	18-20
	G. Recovering plasmids from yeast	20
	H. Obtaining potential positive target plasmids in bacteria	20-21
	I. Determining the specificity of the interaction using mating tests	21-22
<u>VIII.</u>	Appendix	
	A. Polylinker sequences	22
	B. Common false positives	22-23
	C. Plasmid maps	

I. Introduction

The DupLEX-A[™] system is a LexA-based version of the yeast two-hybrid system originally developed by Fields and Song (1). The yeast two-hybrid system has proven to be a powerful tool for identifying proteins from an expression library which can interact with one's protein of interest. The DupLEX-A[™] system was developed as a more versatile and more accurate version of the yeast two-hybrid system (2).

The two-hybrid system of Fields and Song exploits the fact that a yeast transcriptional activator protein, Gal4p, has a separable DNA binding domain and activation domain; neither domain can activate transcription on its own (3). Transcriptional activation is detected only when the binding domain is bound to its DNA recognition sequence and is also tethered to the activation domain. The two-hybrid system involves fusing the Gal4p binding domain with a protein "X" and the Gal4p activation domain with a protein "Y". If "X" and "Y" interact, then a functional Gal4p is restored and transcriptional activation can be detected. If binding sites for Gal4p are placed upstream of a reporter gene (such as lacZ), transcriptional activation can be monitored easily.

The DupLEX-A[™] system utilizes the same basic idea except that the DNA binding protein is the *E. coli* LexA protein while the activation protein is the acid blob domain B42. Neither LexA protein bound upstream of a reporter gene nor B42 alone can activate transcription of the reporter, but if brought together via fusions with two interacting proteins, reporter gene expression can be detected (2).

Advantages of the DupLEX-A[™] system over other yeast two-hybrid systems include:

- A. reduction in the number of false positives obtained since prokaryotic (LexA and B42) rather than eukaryotic (Gal4p) proteins are used
- B. ability to screen potentially toxic target proteins since their expression is galactose-inducible
- C. ability to demonstrate a potential positive's interaction with bait is dependent upon expression of the potential positive
- D. ease of doing a coimmunoprecipitation assay of bait and potential positive since antibodies to HA tag (fused downstream of B42) are available (if antibody to the bait protein is also available)
- E. reporters with varying sensitivities are available so that baits which activate transcription on their own can potentially still be assayed simply by using a less sensitive reporter
- F. ease of determining whether or not a particular bait protein will enter the yeast nucleus and bind LexA operators

References:

1. Fields, S. and Song, O. *Nature* **340**:245-247 (1989).
2. Gyuris, J., Golemis, E.A., Chertkov, H., and Brent, R. *Cell* **75**:791-803 (1993).
3. Chien, C.T. et al. *Proc. Natl. Acad. Sci. USA* **88**:9578-9582 (1991).

II. List of components, additional materials required

NOTE: You will not need all of the components listed below. Read the manual carefully to determine which components will best suit your needs.

A. Yeast strains

EGY48 **MAT** *trp1 his3 ura3 leu2::6 LexAop-LEU2* (high sensitivity)
 EGY194 **MATa** *trp1 his3 ura3 leu2::4 LexAop-LEU2* (medium sensitivity)
 EGY188 **MATa** *trp1 his3 ura3 leu2::2 LexAop-LEU2* (low sensitivity)
 EGY40 **MAT** *trp1 his3 ura3 leu2::0 LexAop-LEU2* (negative control)
 RFY206 **MATa** *trp1 ::hisG his3 200 ura3-52 lys2 201 leu2-3*(mating strain)

Streak on YPD plates and grow at 30°C for 2-3 days; start cultures from single colonies.

B. Reporter gene (lacZ) plasmids (2 ug each)

pSH18-34 *URA3*, 2 um, Ap^R, 8 ops.-lacZ (high sensitivity)
 pJK103 *URA3*, 2 um, Ap^R, 2 ops.-lacZ (medium sensitivity)
 pRB1840 *URA3*, 2 um, Ap^R, 1 op.-lacZ (low sensitivity)
 pJK101 *URA3*, 2 um, Ap^R, GAL1-2 ops.-lacZ (used in repression assay)

C. Bait plasmids

pEG202 *HIS3*, 2 μm, Ap^R (constitutive ADH promoter expresses LexA and is followed by a polylinker for making the bait fusion protein); 10 ug
 pEG202-NLS *HIS3*, 2 μm, Ap^R (similar to pEG202 but with SV40 nuclear localization sequence between LexA and polylinker); 2 ug
 pNLexA *HIS3*, 2 μm, Ap^R (similar to pEG202 except that the LexA sequence is 3' rather than 5' of the polylinker); 2 ug

D. Target plasmid

pJG4-5 *TRP1*, 2 μ m, Ap^R (inducible GAL1 promoter expresses B42-HA tag and is followed by a polylinker for making target fusion protein expression libraries from cDNA); 10 μ g

E. Control plasmids (2 μ g each)

pRHF1 *HIS3*, 2 μ m, Ap^R (ADH promoter expresses LexA- homeodomain of bicoid fusion; used as a positive control in the repression assay and a negative control in the DupLEX-A™ screen)

pSH17-4 *HIS3*, 2 μ m, Ap^R (ADH promoter expresses LexA-GAL4 activation domain; used as a positive control in the DupLEX-A™ screen)

pEG202-Max expresses LexA-Max fusion constitutively; used as a negative control when testing isolated target proteins or as a positive control in the repression assay

pBait constitutively expresses a LexA-bait fusion protein that interacts with the fusion protein from pTarget (see below); can also be used as a negative control when testing isolated target proteins or as a positive control in the repression assay)

pTarget expresses (galactose-dependently) a B42-target fusion protein that interacts with the fusion protein from pBait (see above)

F. Primers (10 μ g each)

[Note: The concentration of each oligo is approximately 80 μ M.]

5' bait fusion primer: 5'-CGT CAG CAG AGC TTC ACC-3'
(used to determine the sequence of the junction between LexA and the bait)

5' target fusion primer: 5'-CTG AGT GGA GAT GCC TCC-3'
(used to determine the reading frame and identity of positive clones; also, can be used with 3' target primer to amplify clone by PCR)

3' target fusion primer: 5'-GCC GAC AAC CTT GAT TG-3'
(used to determine the identity of positive clones; also, can be used with 5' target primer to amplify clone by PCR)

G. Other items included in the kit

pJG4-6 *TRP1*, 2 μ m, Ap^R (similar to pJG4-5 but without B42 activation domain; used to express an isolated target protein in yeast); 10 μ g

sonicated salmon sperm DNA, 5 mg/ml (prepared specially for yeast transformation); 10 mg total [CARRIER DNA]

E. coli strain KC8 *pyrF*, *leuB600*, *trpC*, *hisB463*

DupLEX-A™ System Manual

H. Materials required that are not included in the kit

Note: The specific materials listed below are the ones we have tested in the DupLEX-A™ system. Similar items from other sources may be interchangeable.

1.	<u>Yeast growth media</u>	<u>vendor</u>	<u>catalog #</u>
	peptone	BIO 101	4018-512
	agar	BIO 101	4019-512
	yeast extract	BIO 101	4018-012
	yeast nitrogen base w/o amino acids	BIO 101	4027-012
	raffinose	Sigma	R-0250
	dextrose (glucose)	Fisher	D16-1
	galactose (glucose-free)	Sigma	G-0750
	dropout mix (-his -ura -trp -leu)	BIO101	4540-022
	uracil	Sigma	U-0750
	leucine	Sigma	L-5652
	tryptophan	Sigma	T-0254
	histidine	Fisher	BP382-100
2.	<u>Bacterial growth media</u>	<u>vendor</u>	<u>catalog #</u>
	LB broth	Difco	0402-07-0
	magnesium sulfate	Fisher	BP213-1
	potassium phosphate (monobasic)	Mallinckrodt	7100
	potassium phosphate (dibasic)	Fisher	BP363-500
	sodium citrate	Fisher	BP327-1
	thiamine hydrochloride	Fisher	BP892-100
	ammonium sulfate	Fisher	BP212R-1
	kanamycin	Fisher	BP906-5
	ampicillin	Boehringer Mann.	835 269

3. For yeast transformations:

<u>item</u>	<u>vendor</u>	<u>catalog #</u>
lithium acetate	Fisher	AC19984-2500
polyethylene glycol, 3350	Sigma	P146-3
dimethyl sulfoxide	Fisher	D136-1
hydrochloric acid	Fisher	A144-500
tris base	Fisher	BP152-5
EDTA	Fisher	02793-500

4. For bacterial transformations:

<u>item</u>	<u>vendor</u>	<u>catalog #</u>
glycerol	Sigma	G-5516
ElectroMax DH10B (competent cells)	Gibco/BRL	18290-015
SOC medium	Gibco/BRL	15544-018
electroporator	Gibco/BRL	11613-015

5. For rescuing plasmids from yeast:

<u>item</u>	<u>vendor</u>	<u>catalog #</u>
glass beads, acid-washed	Sigma	G8772
triton-X-100	Fisher	BP151-100
sodium acetate	Fisher	S209-500
ethanol	Aldrich	18,738-0
phenol	Fisher	BP1750I-400
chloroform	Fisher	BP1145-1
isoamyl alcohol	Fisher	BP1150-500
SDS	Sigma	L-4509

6. Filter assay for β -galactosidase and yeast X-gal plates:

<u>item</u>	<u>vendor</u>	<u>catalog #</u>
sodium phosphate (monobasic)	Fisher	BP329-500
sodium phosphate (dibasic)	Mallinckrodt	7914
potassium chloride	Mallinckrodt	6858
-mercaptoethanol	Fisher	03446I-100
nylon membrane filters	MSI	N04SP09025
X-gal	Gold BioTech.	X4281C
N,N-dimethyl formamide	Fisher	BP1160-500

III. Recipes for media and solutions

A. YEAST GROWTH MEDIA

1. YPD (rich medium):

20 g peptone
10 g yeast extract
20 g glucose
one pellet (0.1 g) NaOH [if for plates]
20 g agar [if for plates]

Add 1 liter of distilled water and autoclave for 20 minutes.
For plates, cool to 50°C before pouring.

2. YNB-ura-his-leu-trp (selective medium):

1.7 g yeast nitrogen base w/o amino acids*
5 g ammonium sulfate*
0.6 g -his-ura-trp-leu dropout mix
20 g glucose (or 20 g galactose + 10 g raffinose for gal/raff
media)
20 g agar (if for plates)

Add 1 liter of distilled water and autoclave for 20 minutes.
For plates, cool to 50°C before pouring.

*NOTE: Yeast nitrogen base is sold either with or without ammonium sulfate. The particular one listed on page 6 of this manual does not contain ammonium sulfate, so it must be added separately. If you have a bottle of yeast nitrogen base which instructs you to add 6.7 g/L to make media, then it already contains ammonium sulfate and you should not add any more.

3. other YNB (selective) media:

Add the following amounts of reagents to the YNB-ura-his-leu-trp medium described above (before autoclaving) to make the appropriate medium. Filter sterilize and store at 4°C; microwave briefly if a precipitate forms.

trp = 10 ml of 4 mg/ml stock per liter of medium
(0.04 mg/ml final conc.)
ura = 5 ml of 4 mg/ml stock per liter of medium
(0.02 mg/ml final conc.)
leu = 15 ml of 4 mg/ml stock per liter of medium
(0.06 mg/ml final conc.)
his = 5 ml of 4 mg/ml stock per liter of medium
(0.02 mg/ml final conc.)

For example, to make medium lacking only leucine, add 5 ml of 4 mg/ml uracil, 10 ml of 4 mg/ml tryptophan, and 5 ml of 4 mg/ml histidine to 1 liter of YNB-ura-his-leu-trp medium.

4. Yeast selective X-gal media:

amino acid solution as per 3. above
0.6 g -his-ura-trp-leu dropout mix
1.7 g yeast nitrogen base without amino acids*
5.0 g ammonium sulfate*
20 g glucose (or 20 g galactose + 10 g raffinose for gal media)
20 g agar
900 ml distilled water

Autoclave and cool to 65°C. In a separate bottle, autoclave 7 g of sodium phosphate (dibasic) and 3 g of sodium phosphate (monobasic) in 100 ml of distilled water. Mix the two solutions, add 0.8 ml of 100 mg/ml X-gal (in N,N-dimethyl formamide), and pour plates.

*SEE NOTE ON PAGE 8

B. BACTERIAL GROWTH MEDIA

1. LB medium:

Dissolve 20 g of dry LB Broth (Lennox; from Fisher) in 1 liter of distilled water and autoclave for 15 minutes. For plates, add 15 g of agar per liter of medium before autoclaving and cool to 50°C before pouring.

2. LBA medium (ampicillin selection):

Cool the LB medium above to 50°C and add 2 ml of 50 mg/ml ampicillin (in distilled water, filter-sterilized) per liter of medium. Mix.

3. LBK medium (kanamycin selection):

Cool the LB medium above to 50°C and add 5 ml of 10 mg/ml kanamycin sulfate (in distilled water, filter-sterilized) per liter of medium. Mix.

4. Minimal (-trp) medium:

Prepare the following stocks (autoclaved):

- i. 20 % magnesium sulfate
- ii. 4 mg/ml uracil
- iii. 4 mg/ml histidine
- iv. 4 mg/ml leucine
- v. 20 % glucose

Prepare the following filter-sterilized solutions:

- vi. 50 mg/ml kanamycin sulfate
- vii. 1% thiamine hydrochloride

Autoclave the following two solutions separately:

- viii. 15 g agar in 800 ml distilled water
- ix. 10.5 g potassium phosphate (dibasic)
4.5 g potassium phosphate (monobasic)
1 g ammonium sulfate
0.5 g sodium citrate
160 ml distilled water

Mixing: Cool solutions viii and ix to 50°C, mix, and quickly add 1 ml of solution i, 10 ml of solution ii, 10 ml of solution iii, 10 ml of solution iv, 10 ml of solution v, 1 ml of solution vi, and 0.5 ml of solution vii. Mix well and pour plates immediately.

C. SOLUTIONS

1. 10 X TE

50 ml of 1 M tris (pH 7.5) [0.1 M]
10 ml of 0.5 M EDTA [0.01 M]

Add 440 ml distilled water and autoclave for 20 minutes.

2. 10 X LiOAc

51 g of lithium acetate [1 M]

Bring volume up to 500 ml with distilled water, mix until dissolved, and autoclave for 20 minutes.

3. 50 % PEG-3350

250 g polyethylene glycol-3350

Bring volume up to 500 ml with distilled water, mix until dissolved, and autoclave for 20 minutes.

4. 1 X TE/LiOAc

Right before use, mix 1 part 10 X TE, 1 part 10 X LiOAc, and 8 parts sterile distilled water.

5. 1 X TE/LiOAc/PEG

Right before use, mix 1 part 10 X TE, 1 part 10 X LiOAc, and 8 parts 50 % PEG-3350.

6. Z buffer

16.1 g of sodium phosphate (dibasic) [60 mM]
5.5 g of sodium phosphate (monobasic) [40 mM]
0.75 g of potassium chloride [10 mM]
0.246 g of magnesium sulfate [1 mM]
2.7 ml of β -mercaptoethanol [50 mM]

Dissolve in 1 liter of distilled water. DO NOT AUTOCLAVE!

7. 100 mg/ml X-gal

Dissolve 100 mg of X-gal in 1 ml of N,N-dimethyl formamide and store at -20°C.

8. Plasmid rescue solution

85 ml of distilled water
2 ml of triton-X-100 [2 %]
10 ml of 10 % SDS [1 %]
2 ml of 5 M sodium chloride [0.1 M]
1 ml of 1 M tris (pH 8.0) [0.01 M]
0.2 ml of 0.5 M EDTA [0.001 M]

Mix and store at room temperature.

9. 3 M sodium acetate

40.8 g of sodium acetate trihydrate

Bring volume up to 100 ml with distilled water, mix until dissolved, and autoclave for 20 minutes.

10. 70 % ethanol

Mix 350 ml of absolute ethanol with 150 ml of distilled water and store at -20°C.

11. 10 % glycerol

Mix 100 ml of glycerol with 900 ml of distilled water and autoclave for 20 minutes.

12. 50 % glycerol

Mix 250 ml of glycerol with 250 ml of distilled water and autoclave for 20 minutes.


IV. Library information

All DupLEX-A™ Yeast Two-Hybrid System libraries are made in the B42 activation domain-HA tag expression vector pJG4-5. The cDNA is made by oligo d(T) priming and is cloned unidirectionally between the EcoRI and XhoI sites of pJG4-5 (see Appendix for details). The libraries are provided as ready-to-use plasmid DNA and also as plasmid-containing bacterial cells.

V. Working with yeast

The budding yeast *Saccharomyces cerevisiae* is very amenable to genetic and molecular biological methodologies due to its ability to be transformed by foreign DNA, its highly efficient system of homologous recombination, and its relatively rapid rate of growth. Whereas *E. coli* has a generation time of 30-45 minutes, most yeast strains can double in 90-120 minutes. As with *E. coli*, yeast can be grown on plates or in liquid culture. However, antibiotics which work on *E. coli* do not work on yeast, making good sterile technique mandatory when working with yeast. Finally, the optimum growth temperature for yeast is 28-32°C.

VI. FLOW-CHART FOR DUPLEX-A YEAST TWO-HYBRID SCREEN



VII. DupLEX-A™ Yeast Two-Hybrid System Protocol

A. Constructing the bait

1. Using standard recombinant DNA techniques, subclone your bait protein gene in the correct orientation into the polylinker of pEG202 (see Appendix). Design the bait protein gene subcloning such that it fuses in-frame with LexA. We strongly recommend verifying the sequence of the LexA-bait junction to make sure that a LexA-cDNA fusion protein should be made.

NOTE: We highly recommend testing your bait fusion protein in the assays below before performing a full-scale library screen.

B. Testing the autoactivation potential of the bait (lacZ)

Some bait proteins can activate reporter genes on their own, making a two-hybrid system library screen a waste of time, effort, and resources. However, the DupLEX-A system offers some alternatives if this occurs.

To test for autoactivation by your bait fusion, transform yeast strain EGY48 with the following combinations of vectors:

1. pEG202-Bait + pSH18-34 (test)
2. pSH17-4 + pSH18-34 (strong activation)
3. pRFHM1 + pSH18-34 (no activation)

Small-scale Yeast Transformation Protocol

- a. Grow a 5 ml culture of EGY48 in YPD at 30°C with shaking (overnight). Inoculate by picking a colony off of a streaked plate of EGY48.
- b. Measure the OD₆₀₀ of a 1:10 dilution of the overnight culture. Calculate the OD₆₀₀ of the 5 ml culture and use that to inoculate a 60 ml YPD culture to an OD₆₀₀=0.1. Grow at 30°C with shaking.
- c. When the OD₆₀₀=0.5-0.7 (approximately 4-6 hours after inoculation), pellet the cells by spinning the culture at 1500 x g for 5 minutes. Resuspend in 20 mls of sterile distilled water, spin again, and resuspend the pellet in 0.3 ml of 1 x TE/LiOAc. Put 100 ul into each of three sterile 1.5 ml eppendorf tubes.
- d. Boil the carrier DNA 5 minutes and quickly chill on ice. This is necessary for obtaining a maximum efficiency of transformation; however, it should only be done every third or fourth time the carrier DNA is used.

- e. Add 100 ng of each plasmid DNA and 50 ug of denatured carrier DNA to each tube and mix.
- f. Add 0.3 ml of 1 x TE/LiOAc/PEG, mix by inversion, and put the tubes at 30°C (with or without shaking) for 30 minutes.
- g. Add 70 ul of DMSO (dimethyl sulfoxide) to each, mix by inversion, and put at 42-45°C (without shaking) for 15 minutes.
- h. Spin at 10K rpm in a microcentrifuge for 10 seconds, pour off the supernatant, and resuspend each pellet in 0.5 ml of sterile distilled water.
- i. Spread 50-100 ul of each onto separate YNB(glu)-his-ura plates. Incubate at 30°C for 2-3 days.

Streak 4 colonies from each plate onto another YNB(glu)-his-ura plate. Inoculate at 30°C 1-2 days. Perform a lacZ filter assay or replica to YNB(gal)-his-ura + X-gal plates and grow at 30°C overnight.

Filter Assay

Cut a piece of Whatman 3M paper such that it just fits into a 100 mm petri dish. Put it in an empty dish and add 2 ml of 1 mg/ml X-gal (add 20 ul of X-gal in N,N-dimethyl formamide to 2 ml of Z buffer), making sure the filter is completely wet. Place a similarly-cut nitrocellulose filter on the surface of the plate containing the re-streaked yeast, then pull off and freeze at -70°C for 5 minutes. Take out, thaw, and re-freeze at -70°C. Take out, thaw, and place yeast side up on Whatman filter paper (pre-cut above). Incubate at 30°C for 2 hours.

Expected results: The colonies containing pSH17-4 should turn blue, the colonies containing pRFHM1 should not turn blue, and the colonies containing the pEG202-Bait plasmid may or may not turn blue. If they do not turn blue, then the Bait does not autoactivate reporter gene expression and can be used for screening in the yeast strain EGY48. If the clones containing pEG202-Bait do turn blue in the above assay, then the test should be repeated using pJK103 or pRB1840 instead of pSH18-34. If the Bait does not autoactivate in one of these strains, use that strain to perform the screen. If the Bait still autoactivates even in the least sensitive strain, then you must subclone only portions of the gene encoding your protein into pEG202 and test for a portion that does not autoactivate.

C. Testing the autoactivation potential of the bait (*LEU2*)

Since *LEU2* is the reporter used in the initial screen, it is important not to have a high background of colonies arising due to activation of the *LEU2* gene by the bait alone. Also, for some baits, the *LEU2* reporter in EGY48 is more sensitive than the lacZ reporter on pSH18-34. Therefore, the ability of the bait to autoactivate the *LEU2* reporter should be tested before performing a large screen.

- a. Using a sterile wooden applicator stick, transfer a colony of EGY48 containing the bait plasmid into 0.5 ml of sterile distilled water. Vortex. Dilute 100 ul into 1 ml of sterile distilled water. Vortex. This is Dilution 1. Do three more serial 1:10 dilutions (Dilutions 2-4) such that if Dilution 1 is considered “undiluted”, Dilution 2 = 1:10 diluted, Dilution 3 = 1:100 diluted, and Dilution 4 = 1:1000 diluted.
- b. Plate 100 ul of each of Dilutions 1-4 onto YNB(gal)-his plates and onto YNB(gal)-his-leu plates. Incubate at 30°C for 1-2 days. You should see colonies on the -his plates but not on the -his-leu plates. (Note: galactose plates are used in this experiment since that is the carbon source that will be used during the LEU2 selection step of the large-scale screen.) If you do obtain many colonies on the -his-leu plates, then your bait is autoactivating and you should re-do the assays using your bait in strains EGY194 and EGY188. EGY40 is included as a negative control strain. In addition, three different sensitivity lacZ reporter plasmids are included: pSH18-34> pJK103> pRB1840 (most sensitive>least sensitive).

Once you are convinced that your bait fusion can enter the nucleus and bind to LexA operators without autoactivating either of the two reporter genes, then you are ready to perform a large-scale library screen. Note that for an unknown reason, some baits can autoactivate the reporter genes in a large-scale screen even when they did not autoactivate in small-scale tests. Therefore, although we present you with the protocol for doing a large-scale screen, you may first want to do a “medium-scale” screen, maybe one-fifth the size of a large-scale screen, before performing a large-scale screen.

D. Testing the bait’s ability to enter the nucleus and bind LexA operators

The plasmid pJK101 contains a lacZ reporter gene whose expression is driven by the yeast GAL1 promoter. However, two LexA operators have been placed between the GAL1 promoter and the lacZ gene; LexA fusion proteins will bind to these operators and decrease the level of GAL1-driven lacZ expression.

Repression Assay

- a. Do the following transformations into EGY48 (see protocol under section VI.B.):
 1. pEG202-Bait + pJK101 (test)
 2. pEG202-Max + pJK101 (repression)
 3. pJK101 alone (no repression)

Note: Plate transformations 1 and 2 above onto YNB(glu)-his-ura and transformation 3 onto YNB(glu)-ura plates.

b. Streak 4 colonies from each of plates 1 and 2 above onto YNB(gal)-his-ura + X-gal plates and 4 colonies from plate 3 onto YNB(gal)-ura + X-gal plates. After 12-24 hours at 30°C, it should be evident whether or not the bait fusion can bind to the LexA operators. If some level of repression is observed, then you can proceed with the screen. If no repression is seen, then you should test whether the bait is even being made (you can do a Western blot if antibodies against your bait exist). If you see evidence that the protein is being made in yeast, then try making your bait construct in pEG202-NLS (formerly called pJK202; see Appendix). pEG202-NLS contains a nuclear localization sequence just upstream of the polylinker sequence.

E. Performing a large-scale library transformation

1. Grow a 20 ml overnight culture (at 30°C) of EGY48 (or EGY194 or EGY188) containing the lacZ reporter plasmid (pSH18-34, for example) and your bait plasmid. Grow in YNB(glu)-his-ura medium.
2. The next morning, dilute 100 ul of culture into 0.9 ml of water, mix well, and immediately measure the OD₆₀₀ of the dilution. Multiply by 10 to get the OD₆₀₀ of the undiluted culture. Inoculate 300 ml of YPD medium [in a sterile 2 L flask] with enough of the overnight culture to give an OD₆₀₀ = 0.1. For example, if the dilution of the overnight culture had an OD₆₀₀ = 0.5, then the 20 ml undiluted overnight culture would have an OD₆₀₀ = 0.5 x 10 = 5. Since you want the 300 ml culture to start at OD₆₀₀ = 0.1, the amount of undiluted culture needed would be 0.1/5 x 300 ml = 6 ml. Therefore, you would add 6 ml of undiluted overnight culture to 300 ml of YPD medium and grow at 30°C with vigorous shaking. (Note: Without agitation, yeast cells will settle to the bottom of a flask over time, so you should always swirl the flask before removing any culture to ensure the culture is of uniform density.)
3. Shake the culture vigorously at 30°C until the OD₆₀₀ = 0.5-0.7. This should take 4-5 hours.
4. Harvest the cells by spinning at 1500 x g (3K rpm in a GSA rotor) for 5 minutes (room temperature). Pour off the supernatant.
5. Resuspend in 30 ml of sterile distilled water, transfer to a 50 ml sterile conical tube, and spin at 1500 x g (2.5K rpm in an HL-4 rotor) for 5 minutes (room temperature). Pour off the supernatant.
6. Resuspend the cell pellet in 1.5 ml of 1 x TE/LiOAc. Aliquot 50 ul portions into 30 sterile 1.5 ml eppendorf tubes. (Note: A better transformation efficiency is obtained when doing several small transformations rather than one large transformation.)
7. Denature the carrier DNA by placing the tube in a boiling water bath for five minutes and then quickly placing the tube on ice (do this only if the carrier has not been denatured in the last two or three sets of transformations that have been done). Incubate on ice for five minutes.

8. Add 50 ug (10 ul) of carrier DNA and 1 ug of pJG4-5-based plasmid library DNA to each eppendorf tube. Do not use more than 1 ug of library DNA per tube since multiple plasmids can enter the same yeast cell and give confusing results in later analyses.
9. Add 300 ul of 1 x TE/LiOAc/PEG to each tube and mix by inversion. Incubate (without agitation) at 30°C for 30 minutes.
10. Add 40 ul of DMSO to each tube, mix by inversion, and heat shock by incubating at 42-45°C for 20 minutes.
11. Add 0.6 ml of sterile distilled water to each tube and mix by inversion. Dilute 10 ul of one tube into 990 ul of sterile distilled water, vortex, and plate 100 ul of the dilution onto a 100 mm YNB(glu)-his-ura-trp plate.
Plate all of each tube onto 30 separate 24 cm x 24 cm YNB(glu)-his-ura-trp plates (or 300 ul onto each of 100 150mm YNB(glu)-his-ura-trp plates. Incubate all plates at 30°C for 2-3 days, until colonies appear.
12. Calculate the number of transformants obtained by counting the number of colonies on the 100 mm plate. 100 colonies on the plate corresponds to an efficiency of $10^5/\text{ug}$, so $30 \text{ ug} \times 10^5/\text{ug} = 3 \times 10^6$ total transformants. Similarly, 50 colonies on the 100 mm plate corresponds to $5 \times 10^4/\text{ug}$, or 1.5×10^6 total transformants. A saturating screen of a mammalian library requires at least 2×10^6 transformants.
13. Harvest the transformants as follows:
 - a. Soak a microscope slide in ethanol, then air-dry.
 - b. Pipet 10 ml of sterile distilled water onto each 24 cm x 24 cm plate, scrape off the colonies with the long edge of the microscope slide (taking care to use good sterile technique), and pipet the slurry into a sterile disposable centrifuge tube.
 - c. Centrifuge 5 minutes at 1500 x g (2.5K rpm in an HL-4 rotor) at room temperature. Pour off the supernatant, resuspend the pellet in a total of about 75 ml of sterile distilled water, spin again as above, and pour off the supernatant. Resuspend the pellet in an equal volume of sterile distilled water. Estimate the total volume, add half a volume of sterile 50% glycerol, mix, and freeze 1 ml aliquots at -70°C. Frozen stocks will remain good for at least 1 year when stored at -70°C.

F. Screening for potential positive transformants

1. Titer the number of viable cells
 - a. Thaw an aliquot of the frozen yeast transformants and dilute 1:10 with YNB(gal)-his-ura-trp medium. Shake at 30°C for 4 hours to induce the GAL1 promoter.
 - b. Do serial dilutions of the transformants in YNB(gal)-his-ura-trp and plate 100 ul of each dilution onto separate YNB(gal)-his-ura-trp plates. Incubate at 30°C for 2-3 days, until colonies are visible.

c. Calculate the number of colony-forming units (cfu) per frozen aliquot of yeast transformants.

2. Screen for Leu⁺ colonies

a. To fully screen all the transformants, you should screen about 5-7 times the number of original transformants obtained (see section VI.E.11.). Therefore, if you obtained 3×10^6 transformants originally, you would want to plate about 2×10^7 cfu (from step 1 above) now. Thaw enough of the frozen transformants equal to that number of cfu.

b. The number of cfu per 1 ml frozen aliquot can be used to estimate an OD_{600} of the aliquot since $1 OD_{600} = 2 \times 10^7$ yeast cells. That is, $\# \text{ of cfu} / 2 \times 10^7 = OD_{600}$ of the aliquot. Dilute the appropriate amount of frozen transformants with YNB(gal)-his-ura-trp-leu medium, down to approximately 1×10^7 cells/ml [$OD_{600} = 0.5$]. Shake at 30°C for 4 hours.

c. Plate 100 ul aliquots onto 100 mm YNB(gal)-his-ura-trp-leu plates. You will be screening 1×10^6 cfu per plate (do not exceed this number). Incubate at 30°C for 2 days, until colonies appear.

d. Pick colonies onto a YNB(gal)-his-ura-trp-leu master plate and incubate at 30°C for a few days, until colonies appear. Go back to the plates from step c above and do another master plate for colonies which arise three days after plating and another master plate for colonies which arise four days after plating. Colonies from the first two master plates should definitely be characterized further. Colonies that do not appear until about a week after plating are likely to be artifactual and should not be characterized further, unless they are the only colonies obtained.

3. Test galactose growth dependence and lacZ expression of potential positive transformants

a. Since the expression of the target protein is dependent on galactose, any colonies which can activate *LEU2* and grow on -leu medium in the presence of glucose are false positives and should not be further characterized. Colonies that can grow on -leu medium containing galactose but that cannot grow on -leu medium containing glucose are potentially true positives and should be tested for lacZ expression. Re-streak colonies from the master plates to the following types of plates:

YNB(glu)-his-ura-trp-leu
 YNB(gal)-his-ura-trp-leu
 YNB(glu)-his-ura-trp + X-gal
 YNB(gal)-his-ura-trp + X-gal

b. Incubate at 30°C for 1-2 days, until growth occurs. Potential positive transformants will grow on the -leu (gal) plate but not the -leu (glu) plate, and will turn blue on the X-gal (gal) plate but not on the X-gal (glu) plate.

c. If the number of potential positives is small (<50), then all should be recovered and further characterized. If >50 potential positives are obtained, then you should characterize the first 50 that arise and freeze the rest at -70°C in 20% glycerol.

G. Recovering plasmids from yeast

1. To isolate DNA from the potential positives, grow each one in 2 ml of YNB(glu)-trp medium [or any other -trp medium you have available] at 30°C overnight.
2. Spin down 1.5 ml of each in eppendorf tubes for 10 seconds at maximum speed. Pour off the supernatant, vortex to resuspend the pellet in the residual liquid, and add 200 ul of plasmid rescue solution.
3. Add 100 ul of phenol (tris-sat., pH 8.0) and 100 ul of 24:1 chloroform:isoamyl alcohol. Add about 0.3 g of acid-washed glass beads and vortex vigorously for 2 minutes.
4. Spin at 14 K rpm in a microcentrifuge for 5 minutes at room temperature. Carefully remove 200 ul from the top (aqueous) layer and transfer it to a clean eppendorf tube. Add 20 ul of 3M NaOAc, vortex, and add 440 ul of 95% ethanol. Vortex, spin 20 minutes at 14 K rpm in a microcentrifuge, pipet off and discard the supernatant, wash with 70% ethanol, carefully pipet off and discard the supernatant, vacuum dry the pellet, and resuspend the pellet in 5 ul of sterile distilled water. Use 1 ul to transform *E. coli* KC8 cells to TRP⁺ (see below).

H. Obtaining potential positive target plasmids in bacteria

1. The pellet recovered from the yeast cells in the above procedure consists mostly of yeast RNA and genomic DNA. Very little of it is library plasmid DNA. Also, the trp⁻ *E. coli* strain used (KC8) is not very amenable to transformation. Therefore, electroporation is recommended to recover the library plasmid by transforming into KC8 cells.

Making electrocompetent KC8 cells:


- a. Streak KC8 onto a LBK plate and grow at 37°C overnight.
- b. Pick a single colony and inoculate into 5 ml of LBK medium. Grow at 37°C overnight (with shaking).
- c. Use all 5 ml to inoculate a 500 ml culture in LBK medium. Grow with shaking at 37°C until the OD₆₀₀ = 0.5.
- d. Chill the cells on ice for 30 minutes. Spin at 5 K rpm for 10 minutes at 4°C.

- e. Pour off the supernatant and resuspend the pellet in 300 ml of ice-cold 10% glycerol. Spin again as in d above, pour off the supernatant, and resuspend the pellet in 150 ml of ice-cold 10% glycerol.
- f. Spin again as in d above, pour off the supernatant, resuspend the pellet in 10 ml of ice-cold 10% glycerol, spin again as in d above, carefully pipet off and discard the supernatant, and resuspend the pellet in 2 ml of ice-cold 10% glycerol.
- g. Aliquot into 75 ul portions and store at -70°C.
- h. Electroporate 1 ul of the recovered DNA from step VI.G.4. above into 75 ul of competent KC8 cells using the specifications for your particular electroporator. Plate onto LBA plates and grow at 37°C overnight. Colonies arising at this stage contain either the bait, target, or lacZ reporter plasmid.
- i. To select only those colonies which contain the target plasmid, re-streak (or replica) the colonies from the LBA plate to a minimal (-trp) plate. Grow at 37°C overnight.
- j. Do minipreps on at least 2 colonies from each plate since more than one target plasmid can get into a particular yeast cell. Digest with *EcoRI* + *XhoI* to release the insert and run on a gel to check the size.

I. Determining the specificity of the interaction using mating tests

To test whether the potential positives are specific for the particular bait used, you should test them against other baits with which they should not interact. This is most easily done by transforming your recovered library plasmid back into the yeast strain you used in the screening (EGY48, for example), transforming pSH18-34 + a test bait plasmid into a strain of the opposite mating type, and mating the two strains. **If you put your library plasmid into EGY48, then you should put the reporter plasmid and bait plasmid into RFY206. If you put your library plasmid into EGY194 or EGY188, you should put the reporter plasmid and bait plasmid into EGY40. Alternatively, you can transform the strain you used for screening with two of the plasmids, select the transformed cells, and then transform that strain with the third plasmid.**

Example: Using the small-scale transformation protocol, transform EGY48 with your isolated potential positive target plasmid. Similarly, transform RFY206 with pLexA-Max + pSH18-34. Mate some transformants from each plate by streaking the two in a “+” pattern on a YPD plate.



Incubate the plate at 30°C overnight and replica to a YNB(gal)-his-ura-trp plate the next day. Grow at 30°C for 1-2 days. The only cells that should grow are the ones at the intersection of the two streaks. Perform lacZ filter assays to test for lacZ expression.

pLexA-Max, pBait, and pRHFMI are all bait plasmids that should not interact with your target. You should also test your target against the original bait you used to make sure that it still interacts.

Once you have a positive that passes all the specificity tests, you should transform the DNA into a different *E. coli* strain, do a DNA prep, and sequence a portion of your clone. You can then do a database search to see if you can identify your protein. It is imperative to realize that even a clone that passes all the tests could still be a false positive. For example, in some instances clones have been obtained that appear to interact specifically with a certain bait, but it was already known that the two molecules are located in different parts of the cell; therefore, they are not real interactors. That is why the yeast two-hybrid system should be considered a relatively quick and easy method of obtaining the cloned gene for a protein which may interact with your protein of interest. Once you have the clone, you still need to do more to show that it is biologically relevant. Note that pJG4-6 is included in the kit; it is similar to pJG4-5 except that it does not contain LexA. It can be used to express your target protein alone in yeast.

HAPPY FISHING!

VIII. Appendix

A. Polylinker sequences

1. pEG202 polylinker sequence: (unique sites shown)

5'-CTG GAA TTC CCG GGG ATC CGT CGA CCA TGG CGG CCG CTC
 LexA EcoRI BamHI NcoI NotI

GAG TCG ACC TGC AGC-3'

2. pJG4-5 polylinker sequence:

5'-CCC GAA TTC GGC CGA CTC GAG AAG-3'
EcoRI XhoI

3. pNLexA polylinker sequence:

5'-G AAT TCG CGG CCG CCT CGA GGG ATC CAA TTC ATG AAA GCG-3'

4. pEG202-NLS polylinker sequence:

5'- GTG GAA TTC CCG GGG ATC CGT CGA CCT GCA GCC-3'

B. Common false positives

Although the yeast two-hybrid system has proven to be a powerful tool for detecting protein-protein interactions *in vivo*, all versions of it to date have been plagued with the problem of false positives. Many clones isolated in a library screen will interact with a variety of control baits, suggesting that they are simply artifacts of the screening process and do not represent biologically relevant partners to the particular bait used in the screen. To help researchers better recognize these false positives, Dr. Erica Golemis and Dr. Ilya Serebriiskii of the Fox Chase Cancer Center did a survey of 100 yeast two-hybrid library screens that were performed and what results were obtained. (The results listed here are found on the Web at <http://www.fccc.edu/research/labs/golemis/InteractionTrapInWork.html> and have been reprinted here with the consent of Dr. Golemis.) Of the 100 screens done, 54 appeared to produce biologically relevant interactors, 30 were still in progress, 13 gave no real interactors, and three were aborted due to transactivation of reporter gene expression by the bait or because of instability in the bait fusion protein. Excluding the 30 screens that were still in progress, the results show that 77% of the screens produced real interactors, 19% produced no real interactors, and 4% were aborted due to the reasons mentioned above.

Even for the 77% of the screens that produced real interactors, false positives were still a problem. More than 90% of the screens that produced real interactors also gave false positives. The most common types of false positives obtained were heat shock proteins, ribosomal proteins, ferritin, and ubiquitin. Other false positives reported were cytochrome oxidase, mitochondrial proteins, proteasome subunits, tRNA synthase, collagen-related proteins, zinc finger proteins, vimentin, inorganic pyrophosphatase, PCNA, lamins, elongation factors, and cytoskeletal proteins. Of course, some baits will be involved in real, biologically relevant interactions with one or more of the above listed proteins, so use the list as a guide, not a rule. Remember- even clones that seem to be real interactors should be considered to be putative interactors until a biological relevance for the interaction can be established.

NOTE: FOR RESEARCH PURPOSES ONLY! NOT FOR DIAGNOSTIC OR THERAPEUTIC USAGE.

We will replace, at no cost, any product of ours that does not meet our standard product specifications. No other warranties, expressed or implied, are given with our products. OriGene Technologies, Inc. is not liable for any damages due to the use of this product nor are we liable for the inability to use this product.

PLEASE NOTE THAT THIS KIT IS FOR USE BY THE PURCHASER ONLY AND IS NOT TO BE DISTRIBUTED TO THIRD PARTIES WITHOUT THE WRITTEN CONSENT OF ORIGENE TECHNOLOGIES, INC.

